INSTRUCTION BOOK

Thank you for choosing our Digital Pen Type Thermometer Gamma™ designed without glass and mercury on the basis of the most advanced technologies. The device you are holding now allows you to accurately and comfortably measure your temperature just within 60 seconds.

We are sure that having appraised worthily the high quality and reliability of this device you will become a regular user of the products of Gamma Trade Mark.

Before starting to use this device please study the Instruction Book carefully. The Instruction Book offers all information you need to correctly measure your temperature.

ATTENTION!

- The thermometer is suitable exclusively for measurement of the body temperature (see Measuring Methods)!
- The minimum measurement time until the signalling tone [beep] must be maintained without exception!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Therefore never leave the thermometer unattended to children.
- Protect the thermometer from heating (above 60°C) and striking!
- NEVER boil the thermometer!
- For disposal of the long-life battery, please comply with local regulations for special waste.

1. Features Of The Thermometer

2. Turning On The Thermometer

To turn on the thermometer

1. Press ON/OFF Button over the display window

2. A short signalling tone [beep] signals the thermometer ON.

A display test is performed. Then at an ambient temperature of less than 32°C, an «L» and a flashing «C» appear in the upper right corner of the display. The thermometer is now ready for use.

3. Function Test

Each time it is turned on proper functioning of the thermometer is automatically tested. If a malfunction is detected (measurement inaccuracy), a symbol «ERR» appears on the display, and further measurements become impossible. In this case the thermometer must be replaced.

5. Measuring Methods

In the course of a measurement the current temperature is continuously displayed together with the «C» symbol flashing. The signalling tone (beep-beep-beep 10 times) and «C» no longer flashing means that the measured increase in temperature is now less than 0.1 °C in 16 seconds and that the thermometer displays the measured temperature for reading. Nevertheless note that depending on the measuring method you might have to continue measuring for some minutes more (see Measuring Methods)!

To prolong the life of the battery, it is recommended to turn the thermometer off after use by pressing ON/OFF Button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

Measuring method		Approximate time	Normal temperature
In the mouth (oral) Place the thermometer in one of the two pockets under the tongue, to the left or right of the root of the tongue so that the measuring sensor is in good contact with the tissue. Close your mouth and breath evenly through the nose.		60 seconds	35,7 °C - 37,3 °C
In the anus (rectal) This is the most reliable measuring method that is especially recommended for infants and small children. Carefully insert the tip of the thermometer 2 to 3 cm into the anal aperture.		60 seconds	36,2 °C - 37,7 °C
In the armpit (axillary) From a medical point of view, this method does not provide for accurate readings, therefore it should not be used if precise measurements are required.		5 to 10 minutes*	35,2 °C - 36,7 °C

* keep measuring for 5 more minutes after the signaling tone (depending on the ambient temperature)

6. Storage Of Measured Values

The thermometer automatically stores the result of the last measurement. To see it on the display press and hold the ON/OFF Button down for more than 3

seconds when turning the thermometer on; the maximum temperature stored will be displayed. At the same time, a small «M» for memory will appear in the right corner of the display. After 3 more seconds this value will disappear from the display and the thermometer will switch to the normal measurement mode. Then the stored value is automatically erased.

7. Replacement Of Batteries

When the symbol of upside-down triangle ▼ appears in the lower right corner of the display, the battery is exhausted and needs replacing. Replace the exhausted battery by placing a new one with + at the top.

8. Cleaning

To clean the thermometer wipe it off with a moist towel.

9. Technical Specifications

Digital Pen Type Thermometer Model T-50

Type	Maximum thermometer
Measurement range	32,0°C to 43,9°C Temp < 32,0°C: display «L» for low (too low) Temp > 43,9°C: display «H» for high (too high)
Measurement accuracy	± 0,1°C between 34°C and 42°C at an ambient temperature of 18°C to 28°C
Self-test	Automatic internal check at a test value of 37,0°C; if there is a deviation of >0,1°C, «ERR» (error) is displayed.
Display	Liquid crystal display (LCD) with three digits Smallest unit of display: 0,1°C
Signalling tone	For signalling that the thermometer is ready to use and that the temperature increase is less than 0,1°C/16 seconds.
Memory	For storing the last measured value.
Storage temperature	-10°C to +60°C
Battery	1,5/1,55 V

10. Guarantee

You digital pen type thermometer Model T-50 is covered with a 3-year warranty from the date of purchase. The warranty does not cover the battery and the packing.

Any damage caused by improper handling, accidents, non-observance of the Instruction Book or self-maintained alterations made to the device shall not be covered by the warranty.

The warrantee is only valid upon presentation of the warrantee card which was correctly filled in and sealed or cashier receipt certifying the purchase of the thermometer.

The production date: first three figures of the serial number. The first and the second figures mean the week of production and the third figure stands for the year of production.

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